

Remember GORE-TEX® Suture:

The Perfect Close to Your Vascular Procedures

| | | | | Thread Size: | | | | | | | | |
|---|--------------------------------|--------------------------|-------------------|-----------------------|--|----------------------|----------------------|---|---|-------------|-------------|-------------|
| | | | | CV-8 | CV-7 | CV-6 | CV-5 | CV-4 | CV-3 | CV-2 | | |
| | | | | Approx USP Size: | | | | | | | | |
| | | | | 7-0 | 6-0 | 5-0 | 4-0 | 3-0 | 2-0 | 0 | | |
| | Double-armed Needles | Needle Code | Ethicon Reference | Syneture Reference | Thread Length | Catalogue Number | | | | | | |
| Taper Point 3/8 Circle ⊙ | | TTc-7 | BV130-5 | CV-351, MV-135-5 | 24" – 61 cm | 8K10 | | | | | | |
| | | TTc-9 | BV-1, BV175-8 | CV-1, CVF-1, MV-175-8 | 18" – 46 cm 24" – 61 cm 30" – 76 cm | 8J02 8K02 8M02 | 7J02 7K02 7M02 | 6J02 6K02 6M02 | | | | |
| | | TTc-12 | BV | CV-11, CVF | 24" – 61 cm 30" – 76 cm | | 7M12 | 6M10 | | | | |
| | | TTc-13 | C-1 | CV-11, CVF-11 | 18" – 46 cm 24" – 61 cm 30" – 76 cm 36" – 91 cm | 8M04 | 7J06 7M06 | 6J04 6K04 6M04 6N04 | 5K02 5M02 5N02 | | | |
| | | TT-17 | BB | CV-15, CVF-15 | 36" – 91 cm | | | | 5N06 | 4N04 | 3N02 | |
| | Taper Point 1/2 Circle ⊙ | | THc-9 | N/A | CV-20 | 30" – 76 cm | | | 6M24 | | | |
| | | THc-13 | RB-2 | CV-22, CVF-22 | 18" – 46 cm 24" – 61 cm 30" – 76 cm 36" – 91 cm | | 7J10 7M10 | 6J12 6K12 6M12 | 5K12 5M12 5N12 | | | |
| | | TH-18 | RB-1 | CV-23, CVF-23 | 36" – 91 cm 48" – 122 cm | | | 6N16 | 5N04 | 4N02 | 3N04 | |
| | | TH-22 (with pledgets) | SH-1 | CV-25 | 36" – 91 cm | | | | 3201†† 3202† | 3200†† | | |
| | | TH-26 | SH | V-20 | 36" – 91 cm 48" – 122 cm | | | | 5N20 | 4N06 | 3N08 | 2N02 |
| | | TH-35 | MH | V-26 | 36" – 91 cm 48" – 122 cm | | | | | 4N12 | 3N06 | |
| | | | | | | | | | | | 3U06 | |
| Piercing Point 3/8 Circle Multifaceted tip for tough tissue, tapered body ⊙ | | PT-9 | CC | KV-1 | 18" – 46 cm 24" – 61 cm 30" – 76 cm | 8K08 8M08 | 7J04 7K04 7M04 | 6K06 6M06 | | | | |
| | | PT-13 | CC-1 | KV-11 | 24" – 61 cm 30" – 76 cm 36" – 91 cm | 8M06 | 7K08 7M08 | 6K08 6M08 | 5K08 5N08 | | | |
| | | PT-17 | V-4 | KV-15 | 36" – 91 cm | | | | 5N10 | 4N10 | 3N14 | |
| | Piercing Point 1/2 Circle ⊙ | | PH-13 | N/A | N/A | 30" – 76 cm | | | 6M14 | 5M16 | | |
| | | PH-17 | V-5 | KV-5 | 36" – 91 cm | | | | 5N14 | 4N08 | 3N12 | |
| | | PH-24 | V-7 | KV-7 | 36" – 91 cm 48" – 122 cm | | | | | 4N16 | 3N10 | 2N10 |
| | | PH-30 | N/A | N/A | 42" – 107 cm | | | | | | 3P24 | |
| | | PH-35 | V-34 | KV-34 | 48" – 122 cm | | | | | | 3U20 | |

Double-armed Needles

† With ePTFE pledgets (3) †† With ePTFE pledgets (4)

Catalogue number in bold type indicates 1:1 needle to thread ratio to minimize suture hole bleeding.

To order a three-dozen box, suffix the catalogue number with the letter B. To order a one-dozen box, suffix the catalogue number with the letter A. Some catalogue numbers available in A (one-dozen) boxes only.

Telephone Orders: 800.528.8763, 928.779.2771 (US) • 00800.6334.4673 (EU)

Needle to Thread Ratio

| Thread size: | CV-8 | CV-7 | CV-6 | CV-5 | CV-4 | CV-3 | CV-2 |
|------------------|----------------|----------------|---------------|---------------|---------------|-------------|------|
| Approx USP size: | 7-0 | 6-0 | 5-0 | 4-0 | 3-0 | 2-0 | 0 |
| | | | | | | | |
| TTc-7 (1.6:1) | TTc-9 (1.4:1) | TTc-9 (1:1) | TTc-13 (1:1) | TT-17 (1.2:1) | TT-17 (1:1) | TH-26 (1:1) | |
| TTc-9 (1.9:1) | TTc-12 (1.4:1) | TTc-12 (1:1) | TT-17 (1.2:1) | TH-18 (1.2:1) | TH-18 (1:1) | PH-24 (1:1) | |
| TTc-13 (2.2:1) | TTc-13 (1.6:1) | TTc-13 (1.4:1) | THc-13 (1:1) | TH-22 (1.3:1) | TH-26 (1.2:1) | | |
| PT-9 (1.9:1) | THc-13 (1.6:1) | THc-9 (1:1) | TH-18 (1.2:1) | TH-26 (1.6:1) | TH-35 (1.2:1) | | |
| PT-13 (2.2:1) | PT-9 (1.4:1) | THc-13 (1.4:1) | TH-22 (1.7:1) | TH-35 (1.6:1) | PT-17 (1:1) | | |
| | PT-13 (1.6:1) | TH-18 (1.6:1) | TH-26 (1.7:1) | PT-17 (1.2:1) | PH-17 (1:1) | | |
| | | PT-9 (1:1) | PT-13 (1:1) | PH-17 (1.2:1) | PH-24 (1.2:1) | | |
| | | PT-13 (1.4:1) | PT-17 (1.2:1) | PH-24 (1.6:1) | PH-30 (1.3:1) | | |
| | | PH-13 (1.4:1) | PH-13 (1:1) | PH-35 (1.8:1) | PH-35 (1.5:1) | | |
| | | | PH-17 (1.2:1) | | | | |

Minimize suture hole bleeding with 1:1 needle to thread ratio

Commonly Requested GORE-TEX® Sutures

| Proximal Aorto-bifemoral | | | | | |
|---------------------------|--------------------|---------------------|----------------|------------------------|----------------|
| Thread Size | Thread Length (cm) | Taper Point Needles | Catalog Number | Piercing Point Needles | Catalog Number |
| CV-3 | 91 | TH-26 | 3N08 | PH-24 | 3N10 |
| CV-4 | 91 | TH-26 | 4N06 | PH-24 | 4N16 |
| | 91 | TH-18 | 4N02 | | |
| Distal Aorto-bifemoral | | | | | |
| | 76 | THc-13 | 5M12 | | |
| CV-5 | 76 | TTc-13 | 5M02 | | |
| | 91 | TH-18 | 5N04 | | |
| CV-6 | 76 | TTc-12 | 6M10 | PT-13 | 6M08 |
| Chordae Tendinae | | | | | |
| CV-4 | 91 | TH-22 | 3200* | | |
| | 91/122 | TH-18 | 4N02/4U02 | | |
| CV-5 | 91 | TH-22 | 3201* | | |
| | 91 | TH-18 | 5N04 | | |
| Proximal Femoro-popliteal | | | | | |
| CV-5 | 91 | TTc-13 | 5N02 | PT-13 | 5N08 |
| CV-6 | 76 | TTc-9 | 6M02 | | |
| CV-6 | 76 | TTc-12 | 6M10 | | |
| CV-6 | 76 | TTc-13 | 6M04 | PT-13 | 6M08 |
| Distal Femoro-popliteal | | | | | |
| CV-6 | 76 | TTc-9 | 6M02 | PT-9 | 6M06 |
| CV-7 | 61 | TTc-9 | 7K02 | | |
| | 76 | TTc-9 | 7M02 | PT-9 | 7M04 |
| CV-8 | 61 | TTc-9 | 8K02 | | |

| Carotid Endarterectomy | | | | | |
|------------------------|--------------------|---------------------|----------------|------------------------|----------------|
| Thread Size | Thread Length (cm) | Taper Point Needles | Catalog Number | Piercing Point Needles | Catalog Number |
| CV-6 | 76 | TTc-9 | 6M02 | | |
| | 76 | TTc-13 | 6M04 | | |
| Axillo-Femoral | | | | | |
| CV-5 | 91 | TTc-13 | 5N02 | PT-13 | 5N08 |
| | 76 | TTc-9 | 6M02 | | |
| CV-6 | 76 | TTc-12 | 6M10 | | |
| | 76 | TTc-13 | 6M04 | PT-13 | 6M08 |
| AV Access Graft | | | | | |
| | 61/76 | TTc-9 | 6K02/6M02 | | |
| CV-6 | 61/76 | TTc-12 | 6K10/6M10 | | |
| | 61/76 | TTc-13 | 6K04/6M04 | | |
| CV-7 | 61/76 | TTc-9 | 7K02/7M02 | | |
| AV Fistulas | | | | | |
| CV-7 | 46 | TTc-9 | 7J02 | | |
| CV-8 | 46 | TTc-9 | 8J02 | | |

Catalog number in bold type indicates 1:1 needle to thread ratio to minimize suture hole bleeding
* with ePTFE pledgets (4)

For a full listing of available GORE-TEX® Suture, visit goremedical.com.



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States)

goremedical.com

Products listed may not be available in all markets.
GORE®, GORE-TEX®, and designs are trademarks of W. L. Gore & Associates.
Synature is a trademark of Covidien AG or its affiliate.
© 2008, 2009, 2012 W. L. Gore & Associates, Inc. AM3020-EN3 APRIL 2012